


# B-BRAUN VASOFIX SAFETY PUR IV CATHETER 20G 33mm - sterile



Product Code	23734
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	311790
CND	C0101010201
EAN/UPC	4046963343644

## Description

### BBRAUN I.V. CATHETERS - RADIOPAQUE

#### VASOFIX® SAFETY - sterile

I.V. catheter designed to minimize inadvertent needlesticks, made of PUR.

For single use only, with injection port. Two-way cannula needle, polyurethane, radiopaque, with injection valve for intermittent administration of medicines, complete with stopper closing.

- VASOFIX® SAFETY PUR
- B. Braun code: 4269110S-01
- Gauge: 20 G
- Gauge mm: 33
- Colour: Pink
- Minimum order: box of 50

Latex, PVC, DEHP free.

Flashback technology provides confirmation of proper catheter placement through rapid display of catheter and needle flashbacks. Universal cut with large choice of insertion angles and minimal puncture trauma.

With multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, PL, HU, RO, CZ, SK, BG, GR, HR, RS, SI, EE, LV, LT, SE,FI, NO, RU, TR.