

# SCALP VEIN SET 21G L/L - butterfly type, sterile



Product Code	23768
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2340837
CND	A01010202
EAN/UPC	8032764150563
PARAF	900914553

## Description

STANDARD SCALP VEIN SET - Butterfly type 21G Luer Lock - sterile

Scalp vein sets ensuring minimal trauma during penetration thanks to the triple bevelling. Each element of the needle is treated for painless injection.

The needles are in stainless steel AISI 304 with triple bevelling and totally non-toxic and non-pyrogenic.

The scalp vein set needles are provided with soft flexible grained wings. The tube is made of 30 cm PVC with luer-lock.

- Size: 21G - 0.8x19 mm
- Reference colour: green

- Box of 100 pieces
- Individually packaged in blister
- Box and instructions: GB, FR, IT, DE, ES, PT, GR, RU.