

# SCALP VEIN SET 23G L/L - butterfly type, sterile



Product Code	23770
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2340839
CND	A01010202
EAN/UPC	8032764150570
PARAF	900914565

## Description

STANDARD SCALP VEIN SET - Butterfly type 23G Luer Lock - sterile

Scalp vein sets ensuring minimal trauma during penetration thanks to the triple bevelling. Each element of the needle is treated for painless injection.

The needles are in stainless steel AISI 304 with triple bevelling and totally non-toxic and non-pyrogenic.

The scalp vein set needles are provided with soft flexible grained wings. The tube is made of 30 cm PVC with luer-lock.

- Size: 23G - 0.6x19 mm
- Reference colour: turquoise

- Box of 100 pieces
- Individually packaged in blister
- Box and instructions: GB, FR, IT, DE, ES, PT, GR, RU.