



Product Code	23802
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2369445
CND	A020102010201
EAN/UPC	8034055910979
PARAF	961677933

## Description

Insulin syringe no dead space - assembled needle - 1 ml - 25G - 0.5x16 mm

Available with assembled needle or fixed needle both with no-dead space.

Insulin syringes with fixed needle have the needle directly fixed on the cone and have a double protection (needle cap and plunger cap), assuring the sterility of the device.

Non-toxic and non-pyrogenic, made of transparent polypropylene with graduated insulin scale.  
Needle in stainless steel AISI 304 with triple bevelling.