

# INSULIN SYRINGES NO-DEAD SPACE 30G - 0.5 ml



Product Code	23805
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2326522
CND	A02010601
EAN/UPC	8032764150891
PARAF	915308961

## Description

Insulin syringe no dead space - fixed needle - 0.5 ml - 30G - 0.3x8 mm

Available with assembled needle or fixed needle both with no-dead space.

Insulin syringes with fixed needle have the needle directly fixed on the cone and have a double protection (needle cap and plunger cap), assuring the sterility of the device.

Non-toxic and non-pyrogenic, made of transparent polypropylene with graduated insulin scale.  
Needle in stainless steel AISI 304 with triple bevelling.