

SYRINGES 3 PIECES WITHOUT NEEDLE - 20 ml Eccentric LC



Product Code	23817
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2325134
CND	A020102010202
EAN/UPC	8032764151102
PARAF	915858880

Description

- SYRINGES WITHOUT NEEDLE
- ECCENTRIC LUER CONE
- Capacity: 20 ml

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.