

SYRINGES 3 PIECES WITHOUT NEEDLE - 60 ml CAT



Product Code	23820
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2326507
CND	A02010203
EAN/UPC	8032764150532
PARAF	905324149

Description

- SYRINGES WITHOUT NEEDLE
- CATHETER LUER CONE
- Capacity: 60 ml

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.

High capacity polypropylene siringe.
Sterile, single use, latex free syringes for ear and other purposes.