

# SYRINGES 3 PIECES WITH NEEDLE 21G - 2.5 ml Centric Luer Cone



Product Code	23824
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2325029
CND	A020102010201
EAN/UPC	8032764158026
PARAF	972723973

## Description

- SYRINGES WITH NEEDLE
- CENTRIC LUER CONE
- Capacity: 2.5 ml (extended scale: 3 ml)
- 21G - 0.8 x 38 mm

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Needle in AISI 304 stainless steel triple bevelling.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.