

SYRINGES 3 PIECES WITH NEEDLE 21G - 10 ml Eccentric LC



Product Code	23833
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2358085
CND	A020102010201
EAN/UPC	8032764150600

Description

- SYRINGE WITH NEEDLE
- ECCENTRIC LUER CONE
- Capacity: 10 ml
- 21G - 0.8 x 38 mm

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Needle in AISI 304 stainless steel triple bevelling.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.