

SYRINGES 3 PIECES WITH SAFETY NEEDLE 22G - 2.5 ml Centric Luer Cone



Product Code	23835
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2352591
CND	A02010205
EAN/UPC	8034055936146

Description

- SYRINGES WITH SAFETY NEEDLE
- CENTRIC LUER CONE
- Capacity: 2.5 ml
- 22G - 0.7 x 32 mm

High-quality range of medical devices ensuring an optimal safety.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack.

Needle in AISI 304 stainless steel triple bevelling.

Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.