

# BD CONNECTA THREE WAY STOPCOCKS with extension tube 50 cm - 394951



Product Code	23856
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2371104
CND	A0703
EAN/UPC	8023279238563

## Description

### BD CONNECTA™ PLUS 3 - THREE WAY STOPCOCKS

BD code: 394951  
3-way stopcocks with 50 cm extens. cable

BD Connecta™ Plus 3 are 3-way stopcocks for intravenous therapy and pressure monitoring, with rotating Luer lock.

- Single use, sterile.
- high flow capacity
  - low priming volume (dead space): 0.23 ml
  - full 360° rotation
  - lipid resistant: material allows the use of lipids and chemotherapeutic agents