

BD PLASTIPAK SYRINGES WITHOUT NEEDLE - 50 ml Catheter Cone



Product Code	23877
Unit of sale	box of 60 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	1869867
CND	A020102020102
GTIN	30382903008675
GMDN	47017

Description

BD PLASTIPAK SYRINGES WHITOUT NEEDLES - CATHETER CONE

- Capacity: 50 ml
- Model: Catheter cone
- BD Code: 300867