



Product Code	23887
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2628767
CND	C0101010102
EAN/UPC	8023279238877

Description

BD INSYTE™ IV CATHETERS

BD Instaflash™ needle technology incorporates a notched needle, which clinically demonstrated improved first-stick insertion success, reducing painful hit-and-miss insertions.

BD code: 381234
Gauge: 20G
Gauge mm: 1.09x30 mm
Colour: pink

BD Insyte™ IV catheters are intravascular radiopaque catheters with Luer-Lock cone.

Intended to be inserted into patient's peripheral vascular system for short term use to sample blood, monitor blood pressure or administer fluids.

Needle is internally and externally lubricated with triple sharpening (backcut grinding).

BD Vialon™ material softens up to 70% in the vein, enabling longer dwell times and reducing the chance of mechanical phlebitis and infiltration.