

BD INSYTE AUTOGUARD SHIELDED IV CATHETER 20G x30mm - 381834



Product Code	23901
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1915089
CND	C0101010202
EAN/UPC	8023279239010

Description

BD INSYTE™ AUTOGUARD™ SHIELDED IV CATHETERS

BD INSYTE™ AUTOGUARD™

BD code: 381834

Gauge: 20G

Gauge mm: 1.1x30 mm

Colour: pink

BD Insyte™ Autoguard™ shielded IV catheters are over-the-needle, intravascular (IV) catheters that incorporate advanced safety features including a push button that instantly retracts the needle.

Enhanced Safety

BD push-button needle shielding technology is demonstrated to reduce risk of accidental needlestick injuries.

Insertion success in sight

BD Instaflash™ needle technology incorporates a notched needle, which clinically demonstrated improved first-stick insertion success, reducing painful hit-and-miss insertions.

Reduced complications

BD Vialon™ biomaterial softens up to 70% in the vein, enabling longer dwell times and reducing the chance of mechanical phlebitis. Intended to be inserted into a patient's peripheral vascular system for short term use to sample blood, monitor blood pressure, or administer fluids. Available with or without wings.

Multilingual manual and box in GB, FR, IT, ES, DE, PT, NL, DK, SE, NO, FI, LT, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, RU, TR, Arabic.