



Product Code	23926
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1838349
CND	W0202020101
EAN/UPC	0682607611379
PARAF	977810480

Description

Coagulation testing for professional and self use

Reliable

- quick results for Prothrombin Time (PT) in seconds and International Normalized Ratio (INR)
- large INR range of 0.7-7.0
- individually packed strips with 24 months shelf life
- external optional printer available for quick printing of results
- clinical accuracy: ISO 17593 compliant

Simple operation

- results in about 2 minutes
- requires only a single drop of capillary blood or 15 µL of venous whole blood
- insert strip, apply specimen and read results

User friendly

- handheld portable meter easy to use and ideal for physician offices and home testing
- patients ID, password protected results
- operates with 4 AA batteries or standard AC adaptor

Technical Specifications

Methodology: Optical fluorescence
 Time to results: About 2 minutes
 Memory: 200 results with date/time
 Specimen volume: 15 µL
 Specimen type: Fresh capillary or venous whole blood
 INR measurement range: 0.7-7.0
 HCT range: 25-60 %
 Precision Repeatability and total precision: <10%
 Accuracy (capillary whole blood): $Y=1.076X-0.078$, $R=0.924$
 PC interface: Mini USB
 Power source: 4 AA batteries or AC adaptor
 Battery life: >100 tests
 Meter dimension (LxWxH): 152x72x38 mm (6.0"x2.8"x1.5")
 Display dimension (LxW): 53x29 mm
 Weight excluding batteries: 170 g (6.0 oz)

Standard accessories

- 1 quick reference guide
- 1 AC adaptor
- 4 AA batteries
- 1 carrying case
- 1 user manual (GB, IT) on request DE, ES, PT