



Product Code	23932
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2818170
CND	W0201069080
EAN/UPC	0682607610907
PARAF	977810530

## Description

Provides high quality and reliable results for complete Lipid Panel testing in =2 min.

3 in 1 complete Lipid Panel device tests for Total Cholesterol (CHOL), High Density Lipoprotein (HDL) and Triglycerides (TRIG)

Also includes calculated tests:

Low Density Lipoprotein (LDL), CHOL/HDL ratio and Cardiac Risk Assessment

Reliable

- can test whole blood, serum or plasma
- wide hematocrit range of 0-50%
- small printer available for immediate printing of results

Easy to use

- large LCD display for easy reading
- colour coded test devices with engraved test name for easy identification

Supplied with 2 control devices, 4 batteries, carrying case and manual (GB, IT).

Available FR,ES, DE, PT manual in [www.gimaitaly.com](http://www.gimaitaly.com)

Multilingual box: GB, FR, IT, ES, DE, PT, PL, SE, FI, GR, NL, RO, BG, RU, SA.

## Technical Specifications

Detection principle: Reflectance Photometry

Time to results: = 2 minutes

Memory with date/time: 200 results

Type of test: Quantitative

Tests: CHOL, HDL, TRIG

Calculated tests: LDL, CHOL/HDL ratio, Cardiac Risk Assessment

Specimen volume: 3-1 Lipid Panel Device: 35 µL, CHOL Individual Test Device: 10 µL

Measurement units: mg/dL or mmol/L

Measurement range:

- CHOL: 100-500 mg/dL (2.59-12.93 mmol/L)
- HDL: 15-100 mg/dL (0.39-2.59 mmol/L)
- TRIG: 45-650 mg/dL (0.51-7.34 mmol/L)

Pc interface for data transfer: Mini USB Port

Calibration: Code Chip

Languages on the screen: English abbreviations

Power source: 4 AAA or AC Adaptor (Mini USB, 5V dc, 50 mA)

Battery life: = 1,000 tests

Automatic shut off: 5 minutes

Meter dimensions - weight: 137×79×26 mm (5.4" x 3.1" x 1.0") - 145 g (5.11 oz)  
Display dimensions: 50×50 mm (2.0" x 2.0")