



Product Code	24035
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279240351

Description

XC-2000 CENTRIFUGE

Very silent centrifuge with 6 tube positions, provided with one touch programming display panel and electronic counter to allow precise speed setting from 1000 to 4000 rpm (500 rpm increments).
It is supplied with metal holder for 6 tubes of 10 ml and holder for 6 tubes of 15 ml.
Memory recall of latest operation settings.
Adjustable operating time: 1 - 60 min.
Strong ABS structure.
Automatic shutdown when lid opens.

User manual: GB, IT.

Technical Specifications

Power Supply: 220 V - 50 Hz
Power: 80 W
Speed: Approx 1,000 - 4,000 rpm (increased in increments of 500 rpm's)
Time: 1-60 min
Dimensions: 320x290xh 320 mm
Weight: 6 kg
Capacity: 10 or 15 ml x 6
Max. RCF: 1,790 g