



Product Code	24046
Unit of sale	1 pc
Minimum order	1
Type	IVD - Professional Use
Class	A

RDM (NSIS)	2402832
CND	W020101020101
EAN/UPC	8023279240467
PARAF	977810593

## Description

### GIMA URINE ANALYZER WITH BLUETOOTH

Gima urine analyzer is a high-precision instrument with advanced technologies for clinical inspection of urine. It offers easy and quick operation, exact results and good repeatability.

#### Features:

- high-luminance, white LED, features long life and good stability
- large LCD display
- user-friendly interface
- monitoring the whole testing process, and audible-control prompt
- data transmission function.

#### Application:

- convenient point of care testing for doctor
- chronic diseases, long-term monitoring for nephrotic patients
- health examination
- use as telemedicine terminal.

#### Tested parameters

Glucose (GLU), Bilirubin (BIL), Specific gravity (SG), PH, Ketone body (KET), Occult blood (BLD), Protein (PRO), Urobilinogen (URO), Nitrite (NIT), Leukocytes (LEU), Ascorbic acid (VC).

Can be used only with Urinalysis strips 24076-7.

#### PHMS App

Free App to save data, check report and results.

Multilingual manual and box: GB, FR, IT, ES, SE.

## Technical Specifications

Working mode: One-step

Internal language: Italian and English

Display: 2.4" LCD, resolution: 320x240

Data storage: 500 sample data

Principle: RGB tricolor test theory

Communication interface: Micro USB interface, Bluetooth wireless communication

Repeatability: CV =1%

Stability: CV =1%

Record mode: LCD display, flash ROM data storage

Power supply: Built-in rechargeable lithium battery 3.7 V, 1800 mAh

Host computer: DC 5 V, 1 A

Adapter: AC 100-240 V, 50/60 Hz

Dimension - Weight: 126x73.5x30 mm - 0.18 kg