


# COMBI SCREEN mALB/CREA for detection of microalbuminuria



Product Code	24056
Unit of sale	box of 25 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
--------	-----

RDM (NSIS)	2428580
CND	W0101060204
EAN/UPC	4260033718006

## Description

Combiscreen® mALB/CREA - tube of 25

Visual test for detection of microalbuminuria, based on albumin-to-creatinine ratio that provides accurate information on albumin level in urine, regardless of the urine concentration.

The low detection limit of 22 mg/L enables reliable screening of risk patients.

Result is available in 60 seconds and can be easily read from a table provided with the label or automatically interpreted with Urilyzer® 100 Pro (24049).

Ideal for people with chronic conditions, such as diabetes and high blood pressure.

No 24 h sample collection is necessary.

Designed for healthcare professionals.