

# COMBI SCREEN 9+LEUKO PLUS URINE STRIPS - 9+1 parameters



Product Code	24059
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2428580
CND	W0101060204
EAN/UPC	4260033718129

## Description

### COMBISCREEN® PLUS URINE STRIPS WITH ASCORBIC ACID PROTECTION

CombiScreen® 9+Leuko Plus - 9+1 parameters

Urine strips for visual test or for use with Urilyzer® 100, 100 Pro, 500 Pro and Combiscan 100 (24056\*-8\*, 24062-4 only) \*Not suitable for 500 Pro and Combiscan.

- excellent ascorbic acid protection on glucose and blood test pads offers relevant reduction of false negatives due to interference with ascorbic acid (only CombiScreen® Plus)
- high quality water resistant label made with real colour printing guarantees an exact visual evaluation
- distinct colour changes in the clinically relevant range even at low concentration level
- easy to use: results can be read directly on the strip after about 60 seconds (up to 120 for Leucocytes) and compared with the label on the bottle.

Supplied in alluminium tube in single cardboard box.  
For professional use.

Box all European languages.  
Instructions: GB, FR, IT, ES, PT, DE, NL.  
All other European languages instructions available on [www.gimaitaly.com](http://www.gimaitaly.com).

**WARNING: SHORT EXPIRY DATE FOR TESTS.**  
Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Made in Germany.

## Technical Specifications

COMBISCREEN® PLUS URINE STRIPS - with ascorbic acid protection  
(9+1 PARAMETERS)

- Glucose
- Ascorbic Acid
- Ketones
- Protein
- pH-Value
- Blood
- Nitrite
- Leucocytes
- Bilirubin

