

COMBI SCREEN 10SL PLUS URINE STRIPS - 10 parameters



Product Code	24061
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2428580
CND	W0101060204
EAN/UPC	4260033718013
PARAF	903596346

Description

COMBISCREEN® PLUS URINE STRIPS WITH ASCORBIC ACID PROTECTION

CombiScreen® 10SL Plus - 10 parameters

Urine strips for visual test or for use with Urilyzer® 100, 100 Pro, 500 Pro and Combiscan 100 (24056*-8*, 24062-4 only) *Not suitable for 500 Pro and Combiscan.

- excellent ascorbic acid protection on glucose and blood test pads offers relevant reduction of false negatives due to interference with ascorbic acid (only CombiScreen® Plus)
- high quality water resistant label made with real colour printing guarantees an exact visual evaluation
- distinct colour changes in the clinically relevant range even at low concentration level
- easy to use: results can be read directly on the strip after about 60 seconds (up to 120 for Leucocytes) and compared with the label on the bottle.

Supplied in aluminium tube in single cardboard box.
For professional use.

Box all European languages.

Instructions: GB, FR, IT, ES, PT, DE, NL.

All other European languages instructions available on www.gimaitaly.com.

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Made in Germany.

Technical Specifications

COMBISCREEN® PLUS URINE STRIPS - with ascorbic acid protection
(10 PARAMETERS):

Glucose
Ketones
Protein
pH-Value
Blood
Nitrite
Leucocytes
Spec. Gravity

