



Product Code	24062
Unit of sale	box of 150 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2428580
CND	W0101060204
EAN/UPC	4260033711502
PARAF	903596359

## Description

### COMBISCREEN® CLASSIC URINE STRIPS WITHOUT ASCORBIC ACID PROTECTION

CombiScreen® 11SYS - 11 parameters

Urine strips for visual test or for use with Urilyzer® 100, 100 Pro, 500 Pro and Combiscan 100 (24056\*-8\*, 24062-4 only) \*Not suitable for 500 Pro and Combiscan.

- high quality water resistant label made with real colour printing guarantees an exact visual evaluation
- distinct colour changes in the clinically relevant range even at low concentration level
- easy to use: results can be read directly on the strip after about 60 seconds (up to 120 for Leucocytes) and compared with the label on the bottle.

Supplied in alluminium tube in single cardboard box.  
For professional use.

Box all European languages.  
Instructions: GB, FR, IT, ES, PT, DE, NL.

**WARNING: SHORT EXPIRY DATE FOR TESTS.**  
Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

## Technical Specifications

COMBISCREEN® CLASSIC URINE STRIPS - without ascorbic acid protection  
(11 PARAMETERS):

Glucose  
Ascorbic Acid  
Ketones  
Protein  
pH-Value  
Blood  
Nitrite  
Leucocytes  
Spec. Gravity  
Bilirubin  
Urobilinogen