



Product Code	24073
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2420294
CND	W0101060204
EAN/UPC	4260033718518

Description

10 PARAMETER STRIPS - tube of 100 strips

For professional use.

These strips can be used by professionals for visual tests.
Global check/screening for general health condition.

For use in preliminary screening tests for diabetes, liver diseases, hemolytic diseases, urogenital and kidney disorders and metabolic abnormalities.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, GR, NO, SE, FI, DK, RU, CZ, PL, TR, BG, MT, HU, HR, RO, SK, SL, LV, LT, GA, ET and instructions: GB, FR, IT, ES, PT, DE, NL.

Made in Germany.

Technical Specifications

Parameters:

- Glucose
- Ketones
- Protein
- Ph-Value
- Blood
- Nitrite
- Leucocytes
- Spec. Gravity
- Bilirubin
- Urobilinogen