

# VISUAL URINALYSIS STRIPS 11 PARAMETERS - professional strips



Product Code	24077
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2085574
CND	W0101060204

## Description

11 PARAMETERS STRIPS - tube of 100 strips

PROFESSIONAL USE

Professionals strips for testing and analysis of urine specimen with visual results. Strips can also be examined with Gima urine analyzer (24046).

anti-VC interference ability;

global check/screening for general health condition;

rapid results: from 30 to 120 seconds;

high quality colours chart;

cut-off value measures for each parameter.

Multilingual instructions: GB, FR, IT, ES, PT, SE, HU, GR.

## Technical Specifications

### Parameters

Glucose: Testing and reading time 30 s, Cut-off value 2.8 ~ 5.5 mmol/L

Ascorbic Acid: Testing and reading time 60 s, Cut-off value 0.3 ~ 0.6 mmol/L

Ketone: Testing and reading time 40 s, Cut-off value 0.5 ~ 1.0 mmol/L

Protein: Testing and reading time 60 s, Cut-off value 0.15 ~ 0.3 g/L

pH-Value: Testing and reading time 60 s, Cut-off value 0.5

Blood: Testing and reading time 60 s, Cut-off value 5 - 15 Ery/ $\mu$ L

Nitrite: Testing and reading time 60 s, Cut-off value 13 ~ 22  $\mu$ mol/L

Leucocytes: Testing and reading time 120 s, Cut-off value 5 - 15 Leuko/ $\mu$ L

Specific Gravity: Testing and reading time 45 s, Cut-off value 0.005

Bilirubin: Testing and reading time 30 s, Cut-off value 8.6 - 17  $\mu$ mol/L

Urobilinogen: Testing and reading time 60 s, Cut-off value 3.3 - 16  $\mu$ mol/L