

GIMA GLUCOSE MONITOR mg/dL with lancet device - GB,IT,SE,FI



Product Code	24108
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1832226
CND	W0201060102
EAN/UPC	8023279241082
PARAF	980907036

Description

GIMA Glucose monitor - meter only

Accurate and reliable glucose monitors for self testing and professional use. Extensive clinical studies demonstrate the accuracy of Gima glucose monitoring system, powered by Acon, with fresh capillary, venous, neonatal blood samples.

Comply with ISO15197:2013/EN ISO15197:2015.

Languages box and manual: GB, IT, SE, FI.

Technical Specifications

Technology: Biosensor/Electrochemical, Glucose dehydrogenase (GDH-FAD)

Calibration method: No coding

Result calibration: Plasma-equivalent

Test time: 5 seconds

Sample size: 0.6µL

Sample type: Fresh whole blood (capillary, venous, neonatal)

Test sites: Fingertip, forearm, palm

Hematocrit range: 10-70%

Memory storage: 1,000 results with date and time

Test averaging: 7, 14, 30, 60, 90-day averages

Data transfer: USB

Control solution: 3 levels

Automatic shutoff: 2 minutes after last action

Battery: 2x CR 2032 3.0V coin cell batteries

Battery life: 3,000 measurements

Other functions:

- Automatic detection of insufficient sample
- Sample perturbation detection