

GIMA GLUCOSE MONITOR KIT mg/dL - GB,FR,IT,ES,PT,DE



Product Code	24110
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1832273
CND	W0101050205
EAN/UPC	8023279241105
PARAF	980907048

Description

GIMA GLUCOSE MONITOR mg/dL

Glucose monitor kit
Box and manual languages: GB, FR, IT, ES, PT, DE

Accurate and reliable glucose monitors for self testing and professional use. Extensive clinical studies demonstrate the accuracy of Gima glucose monitoring system, powered by Acon, with fresh capillary, venous, neonatal blood samples. Comply with ISO15197:2013/EN ISO15197:2015.

Kit includes 10 strips, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual.

Technical Specifications

Technology: Biosensor/Electrochemical, Glucose dehydrogenase (GDH-FAD)
Calibration method: No coding
Result calibration: Plasma-equivalent
Test time: 5 seconds
Sample size: 0.6µL
Sample type: Fresh whole blood (capillary, venous, neonatal)
Test sites: Fingertip, forearm, palm
Hematocrit range: 10-70%
Memory storage: 1,000 results with date and time
Test averaging: 7, 14, 30, 60, 90-day averages
Data transfer: USB cable - optional
Control solution: 3 levels
Automatic shutoff: 2 minutes after last action
Battery: 2x CR 2032 3.0V coin cell batteries
Battery life: 3,000 measurements
Other functions: Automatic detection of insufficient sample / 2nd sample application allowed / Sample perturbation detection