



| | |
|---------------|-----------------|
| Product Code | 24148 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | IVD - Self Test |
| Class | |

| | |
|---------|---------------|
| EAN/UPC | 8033937232611 |
|---------|---------------|

Description

MULTICAREIN - English

Supplied with nylon bag 5 cholesterol strips, 5 safety lancets, rapid guide (GB, FR, IT, ES, PT, GR, PL, HU) and GB manual.

Multiparameter diagnostic device for the periodical self-monitoring of risk factors: cholesterol, triglycerides and glucose.

Portable device for easy test procedure and quick results.

Accurate with Amperometric Technology for glucose and Reflectometric Technology for cholesterol and triglycerides.

Switchable mg/dL - mmol/L

Multilingual box: GB, FR, IT, ES, GR, PT, PL, HU, RU and selectable internal software: GB, IT.

Technical Specifications

- Glucose electro-chemical

Measuring range: 10-600 mg/dL - 0.6-33.3 mmol/L

Measuring time: 5"

Blood sample: 0.5 µl

- Cholesterol

Measuring range: 130-400 mg/dL - 3-10.2 mmol/L

Measuring time: 30"-60"

Blood sample: 10 µl

- Triglycerides

Measuring range: 50-500 mg dL - 0.56 - 5.6 mmol/L

Measuring time: 30"-60"

Blood sample: 10 µl

- Memory: 500 test results with date and time

- Display: LCD

- Batteries: 2 lithium batteries Cr 2032 3 Volts (32212) suitable for 1 year or 1,000 measurements

- Size - weight: 97x49x20.5 mm - 65 g (battery included)