



Product Code	24149
Unit of sale	1 kit
Minimum order	1
Type	IVD - Self Test
Class	

EAN/UPC 8033937232697

## Description

### MULTICAREIN starter KIT - GB, IT

Supplied in cardboard box with a nylon bag, 10 glucose, 5 cholesterol, 5 triglycerides strips, 25 safety lancets, user manual (GB, IT), rapid guide and QR code with complete manual in 8 languages (GB, IT, FR, ES, PT, GR, PL, HU).

Multiparameter diagnostic device for the periodical self-monitoring of risk factors: cholesterol, triglycerides and glucose. Portable device for easy test procedure and quick results.

Accurate with Amperometric Technology for glucose and Reflectometric Technology for cholesterol and triglycerides.

Switchable mg/dL - mmol/L

Multilingual box: GB, FR, IT, ES, GR, PT, PL, HU, RU and selectable internal software: GB, IT.

## Technical Specifications

- Glucose electro-chemical  
Measuring range: 10-600 mg/dL - 0.6-33.3 mmol/L  
Measuring time: 5"  
Blood sample: 0.5 µl
- Cholesterol  
Measuring range: 130-400 mg/dL - 3-10.2 mmol/L  
Measuring time: 30"-60"  
Blood sample: 10 µl
- Triglycerides  
Measuring range: 50-500 mg dL - 0.56 - 5.6 mmol/L  
Measuring time: 30"-60"  
Blood sample: 10 µl
- Memory: 500 test results with date and time
- Display: LCD
- Batteries: 2 lithium batteries Cr 2032 3 Volts (32212) suitable for 1 year or 1,000 measurements
- Size - weight: 97x49x20.5 mm - 65 g (battery included)