

# SWANN-MORTON CARBON STEEL STITCH CUTTER

## BLADES - sterile



|               |                |
|---------------|----------------|
| Product Code  | 22420          |
| Unit of sale  | box of 100 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I S            |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2370294       |
| CND        | V0102         |
| EAN/UPC    | 8023279224207 |
| GMDN       | 16224         |

### Description

#### CARBON STEEL STITCH CUTTER BLADE - STERILE

- Carbon steel stitch cutter blade

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.