

SWANN-MORTON SCALPELS WITH STAINLESS STEEL BLADE N. 10 - sterile



Product Code	22430
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2390970
CND	V010102
EAN/UPC	5033955005016
GTIN	05033955005016
GMDN	47569

Description

SWANN-MORTON SCALPELS WITH S/S BLADE

- Disposable scalpels - blade 10 - sterile

Plastic handle scalpels with stainless steel blades, individually packed in sterilized envelope in box of 10 pcs.

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.