



Product Code	24400
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2503557
CND	V901001
EAN/UPC	8023279244007
PARAF	975866548
GMDN	17552

Description

Natural, ecologic method to eliminate completely lice and nits without any harm. Ideal for allergy sufferers and children.

- MANUAL LICE COMB - blister

Ergonomic comb with long metal teeth and a detachable cleaning brush with built-in magnifier.
The long, delicately closed teeth of the comb are designed specifically to detach lice and their nits (tiny eggs).
The round form of the touching point respects the growth direction of the hair without harming its texture.

Multilingual blister: GB, FR, IT, ES, PT, DE, RO, GR.

Detachable cleaning brush with built-in magnifier