



| | |
|---------------|----------------|
| Product Code | 24405 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2263181 |
| CND | V9099 |
| EAN/UPC | 8023279244052 |
| GMDN | 17552 |

Description

Natural, ecologic method to eliminate completely lice and nits without any harm. Ideal for allergy sufferers and children.

- ELECTRIC LICE COMB - GB, FR

This innovative product eradicates and kills lice and nits with a quick increasing electrical discharge of 100-250 Volts from the teeth, imperceptible and harmless to humans.

Auto switch off.

Provided of buzzer indicator and LED light.

Supplied with cleaning brush, 1 AAA battery, multilingual manual in a multilingual blister.

Size: 10.36x7.57x2.66 cm