



Product Code	24454
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2805990
CND	V010401
EAN/UPC	6970934250083
UMDNS	10440
PARAF	987472343

Description

Easy to use, comfortable, safety lancets to meet various capillary blood requirements.

- ROUND SAFETY LANCETS
- Colour: Light blue
- Depth of penetration: 1.8 mm
- Gauge: 26G
- Blood Flow: 5-10 μ l
- Finger stick: YES
- Glucose: YES
- Hematocrit: YES
- Bleeding: Medium

- single-use: the build-in-self-destructing design can prevent reuse of the device

- no accidental injury risk: needle is hidden before and after use to prevent accidental injuries

- sterilization guaranteed: the lancet is Gamma sterilized and the needle tip is protected with a sterility tab until its use allowing a long shelf life

- no blood contamination

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, CZ, SE, RO, SK, NL, DK, FI, EE, HR, HU, LT, LV, NO, SI, BG, RU, GR, Arabic.