



Product Code	24498
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279244984
PARAF	981430352

Description

ALCOHOL TESTER

Compact device to measure alcohol concentration in your breath.
Digital/graphic interface and animation, one-hand operation.
‰ BAC values.
Large display to read exact value of your breath (2 colour backlight).
If the result exceeds alcohol limit, the don't drive logo will appear on display.

Supplied in carton box with 3 blowing tubes Ø 7 mm.

Multilingual box and manual: GB, FR, IT, ES, PT.

Technical Specifications

Range: Low & 0 ~ 1.50 ‰ BAC
Accuracy: ±0.2‰ BAC full range under alert level
Battery: 3 x 1.5V AAA batteries (not included)
Display: Digital display with 0.01‰ BAC resolution
Tolerance compensation: Pre-alert function together with bi-colour backlighting
Breath sampling: about 5 seconds continuous breath sample
Cycle time: about 20 seconds
Size: 107x45x25 mm approx.
Weight: 50 g (without batteries)
Alert level calibrated: 0.5 ‰ BAC
Auto power off after 2 minutes