

SWANN-MORTON STITCH CUTTER WITH RETRACTABLE S/S BLADE - sterile



Product Code	22450
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1420670
CND	V0101
EAN/UPC	8023279224504
GMDN	16224

Description

STITCH CUTTER WITH RETRACTABLE S/S BLADE

- Stitch cutter with retractable blade - sterile

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalledblade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.