



Product Code	24513
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	16427
CND	V9009
EAN/UPC	8023279245134
UMDNS	16380
PARAF	977810668
GMDN	61578

Description

SURGILANCE™ SAFETY AUTOMATIC LANCETS - NEEDLES

Colour: Yellow
 Depth of penetration: 1.0 mm
 Entire needle Gauge: 21G
 Needle tip Gauge: 26G
 Blood Flow: 5-10 µl
 Finger Stick: yes
 Glucose: yes
 Hematocrit: -

Multi purpose capillary blood sampling devices for all testing requirements.
 The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable.
 High-speed delivery and penetration method minimizes patient pain and operator error.
 Latex-free, hypo-allergenic, disposable, sterile.

Multilingual box: GB, FR, IT, ES, DE, PT, NL, SE, DK, RU, GR