



Product Code	24520
Unit of sale	box of 20 pcs
Minimum order	1
Type	IVD - Professional Use
Class	B

RDM (NSIS)	2375791
CND	W0105090103
EAN/UPC	8023279245202
PARAF	987098833

Description

STREP-A - device - box of 20 tests

Provides accurate qualitative results in few minutes, with a simple procedure easy to perform and to interpret. Immediate diagnosis allows fast patient treatment. Reagents are stable for the duration of the shelf life even after opening.

It includes 20 devices (single pack), 2 reagents, 20 tubes and sterile swabs, positive/negative control solution.

Multilingual manual and box (GB, FR, IT, ES, PT, GR, DE, SE, RO, SK)

WARNING: SHORT EXPIRY DATE FOR TESTS. Most tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Technical Specifications

Relative sensitivity: 95.1% (89.6-98.2%)
 Accuracy: 97.1% (95.3-98.4%)
 Relative specificity: 97.8% (95.8-99%)

PROCEDURE

- add 4 drops Reagent A, then drops reagent B
- rotate swab 15 times, wait 1 minute, then remove swab while squeezing tube to remove liquid
- place dropper tip to top of tube
- add 3 drops of specimen
- read results after 5 minutes