



Product Code	24524
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
--------	-----

RDM (NSIS)	2397959
CND	W0105010401
EAN/UPC	6936983123860

## Description

### HELICOBACTER PYLORI SELF TEST - box of 1 test - SELF TEST

Rapid test for qualitative detection of Helicobacter Pylori antigens in human feces specimens. It provides results in 10 minutes.

Packed in single box which includes test cassette, specimen collection tube with extraction buffer, stool collection paper and manual in GB, IT, ES, PT, DE, NL, PL.

Relative sensitivity: 97.6%; relative specificity: 97.9% and accuracy: 97.8%.

**WARNING: SHORT EXPIRY DATE FOR TESTS.** Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 month