


CARDIAC MARKER TEST - 3 parameters



Product Code	24525
--------------	-------

Unit of sale	box of 20 pcs
--------------	---------------

Minimum order	1
---------------	---

Type	IVD - Professional Use
------	------------------------

Class	
-------	--

RDM (NSIS)	1691097
------------	---------

CND	W0102160799
-----	-------------

EAN/UPC	8023279245257
---------	---------------

Description

CARDIAC MARKER 3 parameters: CK-MB, Troponin, Myoglobin - box of 20 tests

For professional use only.

It allows to check cardiac status and stroke events in just 15 minutes.

Ideal to keep under control patients affected by cardiopathies or after heart attack or surgery.

Acute Myocardial Infarction-MI (Heart Attack) is responsible for 1 in 5 deaths in the U.S. After 6 hours of AMI onset, the damage to the cardiac muscle is irreversible and permanent.

This illustrates that early and accurate diagnosis of AMI is essential. ECG is the most widely used method of MI diagnosis.

Unfortunately, as many as 25% of AMI patients have a normal ECG.

GIMA offers an array tests for Biochemical and Protein markers of MI, such as Creatine Kinase Isoenzyme MB (CK-MB), Myoglobin (MB) and Cardiac Troponin I (cTnI) in combination.

It includes 20 droppers, 20 cassettes.

Measurement in plasma, serum or whole blood.

Multilingual manual (GB, FR, IT, ES.) For professional use.

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.