



Product Code	24526
Unit of sale	box of 40 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2498736
CND	W0105090102
EAN/UPC	6936983144704

Description

HELICOBACTER PYLORI TEST - box of 40 tests - PROFESSIONAL

Rapid test for qualitative detection of Helicobacter pylori antibodies in human whole blood, serum or plasma.
Easy to perform and interpret.
Quick results in 10 minutes.

It includes 40 test cassettes, 40 droppers, 1 buffer.

Multilingual box and manual: GB, FR, IT, ES, DE, HU.

Relative sensitivity: 96.8%; relative specificity: 93.0% and accuracy: 94.6%.

WARNING: SHORT EXPIRY DATE FOR TESTS. Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.