



Product Code	24527
--------------	-------

Unit of sale	box of 20 pcs
--------------	---------------

Minimum order	1
---------------	---

Type	IVD - Professional Use
------	------------------------

Class	
-------	--

RDM (NSIS)	1808212
------------	---------

CND	W0103020607
-----	-------------

EAN/UPC	8023279245271
---------	---------------

PARAF	922289487
-------	-----------

Description

CRP TEST - box of 20 tests

Measurement of CRP in whole blood, plasma or serum has two purposes:

1) hs-CRP levels help in assess the risk of cardiovascular and peripheral vascular disease. Higher hs-CRP values have highest risk of cardiovascular disease.
 2) CRP levels help identify and keep track of infections and diseases that cause inflammation (lymphoma, lupus, giant cell arteritis, osteomyelitis, rheumatoid arthritis,).

- use fingerstick whole blood or serum/plasma
- results within 7 minutes
- room temperature storage
- built in performance control and reference band

It includes 20 cassette tests with droppers. For professional use.

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.