



Product Code	24533
Unit of sale	box of 25 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	1519294
CND	W0102160202
EAN/UPC	8023279245332
PARAF	990553125

Description

PSA TEST - box of 25 tests - for professional use

The PSA is a rapid test cassette for semi-quantitative detection of prostate specific antigen, that utilizes a combination of colloidal gold conjugate and anti-PSA antibodies.

- Provides accurate results in 5 minutes.
- sensibility, specificity and accuracy 99%
- detects PSA in whole blood, serum or plasma
- cut-off from 3 ng/ml.

Tests are supplied with droppers, buffer, multilingual instructions (GB, FR, IT, ES).

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.