



Product Code	24537
--------------	-------

Unit of sale	box of 30 pcs
--------------	------------------

Minimum order	1
---------------	---

Type	IVD - Professional Use
------	---------------------------

Class	
-------	--

RDM (NSIS)	1704450
------------	---------

CND	W0105099004
-----	-------------

EAN/UPC	8023279245370
---------	---------------

PARAF	977810694
-------	-----------

Description

INFLUENZA A/B TEST - box of 30 tests (3 packs of 10 pcs.)

For professional use only.

Influenza A/B is a rapid in vitro test to detect type A and B influenza antigens (nucleoprotein) extracted from the nasal specimens. The test is designed to detect nucleoproteins of Flu-A core subtypes, independent of H and N surface proteins (including H1N1, H3N2, H5N1, H5N3, H7N3, H9N2).

A sterile nasal swab allows an easy collection of specimen.

Results in 10-15 minutes.

Includes rapid influenza A/B test panel, swab, test tube, extraction buffer.

WARNING: SHORT EXPIRY DATE FOR TESTS.

Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.