



Product Code	24545
Unit of sale	box of 25 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2368914
CND	W0102160511
EAN/UPC	8023279245455

## Description

MULTI DRUG URINE TEST CUP - 7 drugs - box of 25 tests

A rapid test for simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine. Detects following parameters: (MTD/MDMA) + COC + AMP500 + MET500 + THC + MOP.

- easy to use: no urine handling, just insert the key into the active test and read results in 5 min.
- tightly sealed urine cup: no leaking or spilling of urine specimen
- special design: ensures that urine used in the test, is separated from residual urine, that will remain uncontaminated, to be used for a further test
- key separately packed: only technician can take the test.

Professional use. Manual: GB, FR, IT, ES, PT, DE.

Precision: the results gained = 75% accuracy in  $\pm 25\%$  cut-off level specimen and 100% accuracy in negative

Test panel is tightly fastened to the cup, protecting it from contamination

Max shelf life 22 month

Drug / Cut off NG/ML / Detection time (Drug use over detection time not measured)

MOP: 300 / 5-7 days

COC: 300 / 2-5 days

THC: 50 -

AMP: 500 / 2-5 days

MDMA: 500 / 5-7 days

MTD: 300 / 5-7 days

MET: 500 / 3-5 days