



Product Code	24561
Unit of sale	1 pc
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	Yes
RDM (NSIS)	2450831
CND	W02010701
EAN/UPC	8023279245615

Description

DRUG CUP READER with software and printer

Gima professional drug test reader digitally scans DOA urine cups and gets a complete analysis on touch screen.

It displays and prints the qualitative results only.

- automatic test menu recognition (2D barcode)
- in-process system quality control check
- memory card to store test data
- data management/export archived results

Main benefits:

- automation: it eliminates results interpretation variability/subjectivity
- accurate: performance advantages vs. visual interpretation
- user-friendly: fast results and automatic print
- portable design
- LIS connectivity

Built-in printer and software in GB, FR, IT, ES, PT, DE, GR, HR.

Supplied with print paper roll, power cord and adapter, EU adapter, calibration cup, touch pen and user manual in GB, FR, IT, ES, PT, DE, GR, HR.

Technical Specifications

Mode: photoelectric conversion
 Display: 7" LCD touch screen
 Reading time: 1 - 2.5 minutes
 Detection type: qualitative
 Light source: LED
 Connection: LAN port, 4 USB ports
 Printer: built-in thermal printer
 Power supply: DC 12V, 2A
 Operating altitude: 0-2000 m
 Size: 221x215xh 165 mm
 Weight: < 1.5 kg