



Product Code	24584
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2560245
CND	W0105099099
EAN/UPC	6970277517539

## Description

SARS-CoV-2 ANTIGEN TEST KIT - box of 1

Rapid, membrane-based immunoassay test for the qualitative detection of SARS-CoV-2 Nucleocapsid protein antigens in nasal swab specimen.

For self-testing use.

- results in 15 minutes
- storage condition: 2-30°C
- 24 months shelf life from manufacturing date (average shelf life of Gima stock: 4 to 16 months) Diagnostic sensitivity and diagnostic specificity are 100% (98.9-100%) and 99.5% (98.6-99.9%) respectively.

Kit includes: 1 test cassette, 1 sterile nasal swab, package insert, extraction buffer, biosafety bag and multilingual manual and box in GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, NO, GR.