



Product Code	24586
Unit of sale	box of 20 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2560343
CND	W0105099099
EAN/UPC	6936983154055

Description

SARS-CoV-2 ANTIGEN RAPID TEST KIT - box of 20

SARS-CoV-2 antigen rapid test (nasopharyngeal swab) is a rapid chromatographic immunoassay for qualitative detection of SARS-CoV-2 Nucleocapsid protein antigens in swab specimen.

For professional use only.
Results in 15 minutes.

- storage: 2-30°C
- 24 months shelf life from manufacturing date (average shelf life of Gima stock: 4 to 16 months).
- Relative Sensitivity: 97.6% (95%CI: 93.3%~99.5%)
- Relative Specificity: 99.7% (95%CI: 98.3%~99.9%)
- Accuracy: 99.1% (95%CI: 97.7%~99.8%)

Kit includes: 20 test cassettes, 20 sterile nasopharyngeal swabs, package insert, workstation, extraction buffer, procedure card and multilingual manual in GB, FR, IT, ES, PT, DE, SE, GR.