



Product Code	24588
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

RDM (NSIS)	2122218
CND	W0105099099
EAN/UPC	6970277518062
PARAF	983837509

Description

SALIVA COVID TEST - self test - box of 1

The COVID-19 Antigen Rapid Test (Oral Fluid) is a qualitative rapid chromatographic immunoassay for the qualitative detection of SARSCoV-2 Nucleocapsid (N) protein antigens in oral fluid specimens from individuals with suspected infection. It's a single-use test kit intended to detect the coronavirus SARS-CoV-2 through human oral fluid.

This test is designed for self-testing.

Research shows that the SARS-CoV-2 virus is stable in oral fluid for long periods of time at warm temperatures.

- fast results: 15 minutes
 - high accuracy
 - convenient procedure
 - easy interpretation
 - simple operation
 - relative specificity: 99,3%
- Storage temperature: 2-30°C

Shelf life: 12 months (Gima stock may vary from 4 to 10 months)

Standard accessories

Test devices
Collection devices (funnel, tube and tube tip)
Buffer
Package insert
Biosafety bag
Multilingual box: GB, FR, IT, ES, DE, PT.