

DIAMANTINE DISPOSABLE SCALPELS WITH S/S BLADE N. 23 - sterile



| | |
|---------------|----------------|
| Product Code | 22467 |
| Unit of sale | box of 10 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|------------|---------------|
| UK-REP | Yes |
| CH-REP | No |
| RDM (NSIS) | 2391016 |
| CND | V010301 |
| EAN/UPC | 8023279224672 |

Description

DIAMANTINE™ EDGE HIGH PERFORMANCE SCALPELS - sterile

- Scalpels - blade 23
- Sterile

High-Performance Edge blade for extreme and durable sharpness.

- sensory grip handle enhances cutting action and control of the line of incision
- each handle has a graduated 6 cm ruler to help determine incision length (indicative only)
- protective sheath covers blade to help prevent accidental cuts
- latex and phthalates free

Sensory grip handle.
Stainless steel Diamantine Edge high performance coated blade.
Hardness over 2,500 Vickers

Size is clearly marked on blade and handle for easy identification.