

FINGER GONIOMETER 6" - stainless steel



Product Code	24900
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279249002

Description

FINGER GONIOMETER 6" - 15 cm

Finger goniometer measures range of motion of finger joints.
Linear scale in inches and centimeters on arms.
180° scales marked in 5° increments.

Multilanguage manual: GB, FR, IT, ES.

Approximate measurement of values.